

NIST Laboratory Environment Monitoring and Alert System (LEMAS)

Documentation for building, usage, and configuration

Michael Braine, michael.braine@nist.gov

June 11, 2019 (Original: August 22, 2018)

Revision 1.3

What is this, and why should I read it?

This documentation outlines device setup and customization options should a user desire to perform changes to the controls that dictate when a MMS message or email is sent to lab users. It assumes no understanding of Python3 syntax, Python libraries that are used, or Linux commands; however, the reader should be, at a minimum, comfortable on a command line interface and know what a raspberry (raz,berē,'raz,b(ə)rē – an edible soft fruit related to the blackberry, consisting of a cluster of reddish-pink drupelets) looks like. The curator, in the context within, is the current maintainer of all the monitoring devices (primary curator) and could include additional members whom are designated to receive the same messages as the primary curator.

Table of Contents

Components required for a new device

Contents of LEMASdist package

Logic for when a message is sent

Installation of LEMAS with LEMASRaspbian image

Manual installation of LEMAS

Device customization

Components required for a new device

Raspberry Pi 3

Raspberry Pi 3 power cable

microSD card, >= 16 GB preferred

Ethernet cable (optional, alternatives are communications over WiFi or no communications at all)

USB to microUSB cable (optional, power for Adafruit display)

HDMI male to HDMI male cable (optional, video feed for display)

Sensor(s)

Interface cable for sensor(s)

<picture of setup>

Using a SDXC microSD card (typically larger than 16 GB capacity, but includes some 16 GB cards) will require exFAT32 formatting.

The current revision of LEMASrun.py was tested to work with the above components, and this document written with the assumption that the device will be using the above parts list. A Raspberry Pi 3 Model B (hereby RPi for the rest of the text) is not required, but knowledge of Python3 (or Python2, if you prefer that flavor) and the native OS will be required for setup, especially if running the program on boot is desired.

Contents of LEMASdist package

<i>Filename</i>	<i>Description</i>	<i>Raspbian Install Location</i>
LEMASRun.py	Python3 script for logging and alerts	/home/pi/LEMASdist
LEMASRunQuiet.py	Python3 script for logging and no alerts	/home/pi/LEMASdist
LEMASRunLoud.py	Python3 script for logging and alerts	/home/pi/LEMASdist
InstrInterface.py	Sensor communication functions	/home/pi/LEMASdist
ServerInfo.py	Setting up sending messages through SMTP	/home/pi/LEMASdist
messages.py	Custom messages to send under various conditions	/home/pi/LEMASdist
NoContact.list	list of names to not send messages to	/home/pi/LEMASdist
LabSettings.py	user-configurable settings file	/home/pi/LEMASdist
LabID.py	identifies the lab that the device resides	/home/pi/LEMASdist
SensorSerial.py	serial number for sensor pair with device	/home/pi/LEMASdist
corrections.py	list of corrections for all sensors	/home/pi/LEMASdist
testmsgdate.csv	date to send test message to users	/home/pi/LEMASdist
Contacts.py	list of contact information and lab users	/home/pi/LEMASdist
Tcontrols.py	temperature controls for each lab	/home/pi/LEMASdist
Rhcontrols.py	humidity controls for each lab	/home/pi/LEMASdist
LEMAS	shell commands for starting, making quiet/loud	/home/pi/.bin
setup-LEMAS	shell command to configure files for distribution	/home/pi/LEMASdist
install-LEMAS	shell command to install configured files to RPi	/home/pi/LEMASdist
launch-LEMAS	executable Bash script to launch system	/home/pi/Desktop
config.txt	configuration for Adafruit display	/boot
lightdm.conf	configuration for disabling screensaver	/etc/lightdm
README	github readme	/home/pi/LEMASdist
LICENSE	generic NIST software license	/home/pi/LEMASdist
version	details of the installed LEMAS version	/home/pi/LEMASdist

Logic for when a message is sent

A message is sent from any single device under one of three conditions:

- 1) When the set environment has exceeded limits, with a graph
- 2) When the set environment has returned within limits for a period of time, with a graph
- 3) When the set environment has exceeded limits, and the environment (T or RH) changes in increments of TincSet or RHincSet (LabSettings.py) until the environment returns within limits for a period of time, with a graph

In the first case, the device will send a message with an attached graph to its laboratory owners should temperature or humidity limits be exceeded. These limits can be easily changed, detailed in section *Configuring the device for a lab with existing settings*.

In the second case, the device will send a message with an attached graph to its laboratory owners should temperature or humidity conditions return within limits after a period of time. These limits are the same limits that decide when a message is sent under the third condition and can be easily changed, detailed in section *Configuring the device for a lab with existing settings*.

In the third case, when the environment is outside of the setpoints, the device will send a message for each occasion the environment changes by TincSet or RhincSet. E.g., the temperature upper limit is 25 deg. C, and TincSet is 5 deg. C, there will be a message generated at 30 deg. C, 35, deg. C, etc, until the temperature returns below 25 deg. C.

Installation of LEMAS with LEMASRaspbian image (NIST only, unless hosting service is found)

The LEMASRaspbian image is a copy of the entire ext4 filesystem of a working LEMAS device. The Raspbian kernel is configured and compiled to be securely used on the NIST network. LEMAS is already installed in the image.

1. With a microSD card writer, copy the latest image onto a blank microSD card.
2. Upon supplying power to the Raspberry Pi, the OS will load and LEMAS will launch by itself.

Manual installation of LEMAS

Device setup from a blank microSD card

Required to install Raspbian to a microSD card is a microSD card reader, or adapter for SD-to-microSD.

Installing Raspbian OS from a Windows Machine

1. Navigate a web browser to <https://www.raspberrypi.org/downloads/raspbian/> and download the latest full Raspbian ZIP image for Raspberry Pi 3.
2. Navigate a web browser to <https://etcher.io>, download and install the latest Etcher release.
3. Insert the microSD card into the computer. Launch etcher. Select the downloaded Raspbian image file, select the drive containing the microSD card, select Flash!.
4. Once Raspbian has been written to the microSD card and verified by Etcher, remove the microSD card from the computer and insert it into the Raspberry Pi. Connect the Raspberry Pi to a computer monitor (*NOT THE ADAFRUIT DISPLAY*, the RPi must be configured because the display does not have onboard video codec feedback. Running the installer BASH will configure the Adafruit display). Supply power to the Raspberry Pi to through the micro USB cable. The RPi will boot directly into the default user “pi” with default password “raspberrry”.

Updating Raspbian OS and Python3

1. Activate the US-based keyboard layout by navigating through the start menu (raspberry icon > Preferences > Mouse and Keyboard). Click the Keyboard tab, click on the Keyboard Layout button, select United States under Country, and English (US) under Variant.
2. Connect the device to the internet with Ethernet and use the network icon on the taskbar to connect to the NIST Research Network. Ask the OU IT Security Officer or group IT Security Representative for assistance if needed. Do note that you will need access to repositories external to NIST at least once when setting up the REN connection. Alternatively, the NIST Guest WiFi is suitable for updates, since the Raspberry Pi 3 Model B has a built-in WiFi chip. Similarly, use the network icon to connect to NIST Guest WiFi.
3. In Raspbian, open a terminal either by clicking the button on the taskbar or navigating through the start menu (raspberry icon > Accessories > Terminal). Type the following command in the terminal and press enter:
`sudo apt-get update`
 This command will ask for a list of updates of all installed packages that are managed by the apt-get package manager, including Raspbian OS and Python 3. Once the list has been updated, type the following command in the terminal and press enter:
`sudo apt-get upgrade`
 This command will download and install the updates. Reply y to any y/n prompts.
4. (optional, for freeing extra disk space) In the terminal, enter:
`sudo apt-get remove minecraft-pi -y && sudo apt-get remove bluej -y && sudo apt-get remove geany -y && sudo apt-get remove greenfoot -y && sudo apt-get remove nodered -y && sudo apt-get remove scratch -y && sudo apt-get remove scratch2 -y && sudo apt-get remove sense-hat -y && sudo apt-get remove sonic-pi -y`
 These commands will remove unused programs preloaded on Raspbian.
5. In the terminal, enter:
`sudo apt-get autoremove`
`sudo apt-get autoclean`
 These commands will clean up any files that may no longer be used due to the update processes.
6. In the terminal, enter:
`sudo apt-get install python3-pyserial python3-numpy python3-matplotlib`
 This command installs the PySerial, NumPy, and matplotlib libraries for Python3 for serial interfacing, numerical math operations similar to Matlab, and matrix plotting library, respectively. Note that the python3 designation here is not a Linux command to run Python, it is asking for the libraries in the repositories to download that are built for Python3.
7. In the terminal, enter:
`sudo python3 /usr/bin/pip install minimalmodbus`
 This will install the minimalmodbus library for Modbus RTU communications with Python's pip installer for Python 3.

Configuring the Raspberry Pi 3

8. On a flash drive, load the LEMASdist folder provided by the curator. Plug the flash drive into the RPi on any USB slot. Open file manager with the folder icon on the taskbar, or navigating through the start menu (raspberry icon > Accessories > File Manager).
9. If the dropdown menu does not already read "Places", click the dropdown that reads "Directory Tree" and change to "Places". Access the flash drive and copy/paste the LEMASdist folder onto the desktop.
10. In the terminal, use the copy command (cp) to move the config.txt file to /boot/:
`sudo cp home/pi/Desktop/LEMASdist/config.txt /boot/config.txt`
 The config.txt file will allow the Adafruit display to work with the RPi upon restarting with the display connected and powered. If tilde (~) does not appear, rather another symbol, then the keyboard needs changed to English (US) (see Step 9).
11. In the terminal, enter the following command to change audio output from over HDMI to the 3.5mm jack. This will prevent potential issues when using the Adafruit display:
`sudo amixer cset numid=3 1`

12. In the terminal, enter the following command:

```
sudo cp home/pi/Desktop/LEMASdist/lightdm.conf etc/lightdm/lightdm.conf
```

This copies a file which will deactivate the blank screen screensaver so that graphed results can always be viewed without mouse/keyboard action.

13. Add a lxterminal session to launch at user login by opening the autostart file. In a terminal:

```
sudo echo '@lxterminal -e /usr/bin/python3 /home/pi/LEMASdist/LEMASRun.py' >> /etc/xdg/lxsession/LXDE-pi/autostart
```

This will instruct the system to open a terminal and execute the logging program upon user login.

The Raspberry Pi will begin logging and sending messages once on the REN upon restart; however, the Python logging script will need configured for the lab it will be monitoring in order for the device to identify itself correctly and send messages under the desired outage conditions. See sections *Configuring the device for a lab with existing settings*, *Adding new labs and new users*, and *Adding new lab settings*, as appropriate.

Setting file ownership and permissions

Open a terminal (CTRL+ALT+T) and use the following commands to change file ownership (chown) and file permissions (chmod):

```
sudo chown -R pi:pi /home/pi/LEMASdist/
```

```
sudo chown -R pi:pi /home/pi/Desktop/
```

```
sudo chmod u=rwx,g=rx,o=rx /home/pi/Desktop/launch-LEMAS
```

For NIST

If the device data will populate a website, the site maintainer must have an account on the pi with permission to access the data. Since the script is run locally by user pi on startup, the account must be given permission to access group pi.

Typically, the account is called dmgenenv:

```
sudo useradd dmgenenv
```

Give the login credentials to the site maintainer. Then, add user dmgenenv to group pi:

```
sudo usermod -a -G pi dmgenenv
```

user dmgenenv can now access data generated by the scripts run by user pi.

Device Customization

Configuring the device for a lab with existing settings

If a configuration does not exist for the lab, then the following steps must be skipped and section *Adding new lab settings* followed first.

1. In the terminal, execute the following command:

```
leafpad ~/LEMASdist/LabID.py
```

This command opens the file containing lab settings in the leafpad text editor.

2. In the leafpad text editor, navigate to the line of text that reads: labID = '219/G032'. This line will need updated to reflect the lab that the device will be monitoring. The format of the labID variable is, inside single quotes ('') <building number>/<lab number>. For example, if the device will be used in building 220, room B114, the line must read:

```
labID = '220/B114'
```

Once the change is made, press Ctrl+S to save the change, then Ctrl+Q to quit leafpad.

3. Properly reboot the RPi by performing Start > Restart, or opening a terminal (CTRL+ALT+T) and using the command `sudo reboot`.

The Raspberry Pi will now begin logging upon boot.

Changing existing lab controls

If a configuration exists for a lab, one must only edit `Tcontrols.py` or `RHcontrols.py` for either temperature controls or humidity controls, respectively. Both have an identical procedure. For example, temperature controls:

1. Open `Tcontrols.py` in a text editor:

```
leafpad ~/LEMASdist/Tcontrols.py
```
2. Find the `Tcontrols` dictionary that contains the desired lab. Within the brackets, modify the values to the desired setting. First is minimum, second is maximum. E.g. for a device with labID 220/B114:

```
Tcontrols['220/B114'] = [19.5, 20.5]
```

Adding users to an existing lab

If a configuration exists for a lab, one must only edit `Contacts.py`. Note that `Contacts.py.gpg` must first have been decrypted into `Contacts.py`.

1. Open `Contacts.py` in a text editor:

```
leafpad ~/LEMASdist/Contacts.py
```
2. Find the `labusers` dictionary that contains the desired lab. Add a comma after the last user, before the closing square bracket (`]`) and add the name of the user in single quotes. Their name spelling and capital letter usage must match their name in the `allcontacts` dictionary.

Adding new labs and assigning contacts

If the contacts already exist, then section the *Adding New Contacts* can be ignored. Otherwise, the contacts need added and *Adding new contacts* must be followed as well. Note: `Contacts.py.gpg` must first have been decrypted into `Contacts.py`.

1. If `Contacts.py` is not already opened in leafpad, in a terminal, execute the following command:

```
leafpad ~/LEMASdist/Contacts.py
```

This command opens the lab monitoring program in the leafpad text
2. To add a new lab, navigate to the line that reads:

```
#labusers['<building>/<room>'] = ['user1', 'user2']
```

This is the line that needs edited. Delete the symbol for comment line (`#`) and replace the text inside the brackets (`<>`), including the brackets, with the necessary information. The names entered for users must match those in the `allcontacts` dictionary.

Adding new contacts

If a lab does not exist for the contacts, follow section *Adding new labs* first. Note: `Contacts.py.gpg` must first have been decrypted into `Contacts.py`.

1. Open `Contacts.py` in leafpad:

```
leafpad ~/LEMASdist/Contacts.py
```

2. To add a new contact, navigate to the line that reads:
`#allcontacts['user1'] = ['<number>@<carrier>', '<email>@<domain>']`

This is the line that needs edited. Delete the symbol for comment line (#) and replace the text inside the brackets (<>), including the brackets, with the necessary information. The spelling and capital letter usage for the user must match that in the labusers dictionary.

Adding new lab settings

This section is intended for when a device is added to a lab for which there are no environment settings. The settings are contained in Tcontrols.py and Rhcontrols.py, set as the Tcontrols and Rhcontrols variables, respectively. Adding users to a lab, or adding a new lab, is detailed in section *Adding new labs and new users*. The new lab should be added first.

1. In the terminal, execute the following command:
`leafpad ~/LEMASdist/Tcontrols.py`
to edit temperature controls or:
`leafpad ~/LEMASdist/Rhcontrols.py`
edit humidity controls.
2. On the last line of the file, enter a new line below. Enter into the line the format for Tcontrols:
`['<building>/<room>'] = [<minimum temperature>, <maximum temperature>]`
3. Remove any less, than, greater than symbols. For example, for a lab in building 219, G032 with minimum temperature of 19 deg. C and maximum temperature of 21 deg. C, the line will read as follows:
`Tcontrols['219/G032'] = [19, 21]`
4. On the last line of Rhcontrols.py, enter a new line below. Enter into the line the format for RHcontrols
`Rhcontrols['<building>/<room>'] = [<minimum humidity>, <maximum humidity>]`
5. Remove any less, than, greater than symbols. For example, for a lab in building 219, G032 with minimum humidity of 30 %RH and maximum humidity of 60 %RH, the line will read as follows:
`Rhcontrols['219/G032'] = [30, 60]`

Setting a date and time for a test message

This section is intended for assigning a date and time for sending a test message to all users. This assumes the intended device is set up with users and outage conditions.

1. In Raspbian, on the device to send a test message from, open a terminal. Use the following command to open the configuration file:
`sudo leafpad /home/pi/LEMASdist/testmsgdate.py`
2. Edit the TestmsgDate variable with the date and time, in 24-hour format, a test message is desired to be sent. Surround the date and time with single quotes (''). Use the following examples to input the date and time in the correct format:

`TestmsgDate = 'February, 01 2118 18:30:00'`

or

`TestmsgDate = 'August, 22 1986 06:30:00'`

Both of these settings will send a test message to all users of the device at 6:30 PM or 6:30 AM, respectively, on their respective month, day, and year. The latter will not send any messages at all because the date has already passed...hopefully. If it hasn't, then a time paradox must have occurred because the author of this document and logging program has not yet been born. It is very important to note that setting the date and time as February 1, 2118 18:30:00 may not have the intended results. Notice the absent 0 in the day of the month. Do not enter any other lines, there must only be the first line containing the date and time. This format *must* be followed in order for the date and time to be parsed correctly. For reference, the strftime format is: "%B %d, %Y %H:%M:%S"

The website: <http://strftime.org/> details the meaning of each symbol.

Appendices

Appendix A: Common carrier gateways

See <https://martinfitzpatrick.name/list-of-email-to-sms-gateways/> for a more complete list.

# Carrier	Email to SMS Gateway
# Alltel	[10-digit phone number]@message.alltel.com
#	Example: 1234567890@message.alltel.com
# AT&T (formerly Cingular)	[10-digit phone number]@txt.att.net
#	[10-digit phone number]@mms.att.net (MMS)
#	[10-digit phone number]@cingularme.com
#	Example: 1234567890@txt.att.net
# Boost Mobile	[10-digit phone number]@myboostmobile.com
#	Example: 1234567890@myboostmobile.com
# Nextel (now Sprint Nextel)	[10-digit telephone number]@messaging.nextel.com
#	Example: 1234567890@messaging.nextel.com
# Sprint PCS (now Sprint Nextel)	[10-digit phone number]@messaging.sprintpcs.com
#	[10-digit phone number]@pm.sprint.com (MMS)
#	Example: 1234567890@messaging.sprintpcs.com
# T-Mobile	[10-digit phone number]@tmomail.net
#	Example: 1234567890@tmomail.net
# US Cellular	[10-digit phone number]@email.uscc.net (SMS)
#	[10-digit phone number]@mms.uscc.net (MMS)
#	Example: 1234567890@email.uscc.net
# Verizon	[10-digit phone number]@vtext.com
#	[10-digit phone number]@vzwpx.com (MMS)
#	Example: 1234567890@vtext.com
# Virgin Mobile USA	[10-digit phone number]@vmobl.com
#	Example: 1234567890@vmobl.com

Appendix B: Locally compiling custom Raspbian Image (NIST-only)

The Raspbian kernel must have additional features enabled to be allowed on the NIST network. For reference, the instructions for custom building a Raspbian kernel is found on the Raspberry Pi website (<https://www.raspberrypi.org/documentation/linux/kernel/building.md>). The kernel should be built on a Raspberry Pi 3 for the reason that weird behaviors were noted to occur when cross-compiled with a Ubuntu 14 system. An internet connection is required to grab the OS source code and build dependencies.

Open a terminal with CTRL+ALT+T. Install git and build dependencies using:
sudo apt install git bc

In the terminal, obtain the source code using:
git clone -depth=1 <https://github.com/raspberrypi/linux>

Compute the default build configuration:

```
cd linux
```

```
KERNEL=kernel7
```

```
make bcm2709_defconfig
```

Use make-menuconfig to configure the kernel to enable CONFIG_CPU_SW_DOMAIN and CONFIG_AUDIT. Enter make-menuconfig in the terminal to launch menuconfig, and search for the required features in the kernel menu. Save when finished.

Build and install the kernel using:

```
make -j4 zImage modules dtbs
```

```
sudo make modules_install
```

```
sudo cp arch/arm/boot/dts/*.dtb /boot/
```

```
sudo cp arch/arm/boot/dts/overlays/*.dtb* /boot/overlays
```

```
sudo cp arch/arm/boot/zImage /boot/$KERNEL.img
```

Upon rebooting, the Raspberry Pi will use the new, configured kernel.